



A general purpose thermal roll printer, the ATP2 is ideal for use with most balances and scales. Designed to be efficient and compact, the printer occupies minimal space on a bench or tabletop. Operation is simple, with an on/off switch, paper feed button and two LED indicators. Ideal for printing applications, the ATP2 is able to perform a variety of other commands, such as setting tabs, spacing and carriage return. Users can specify absolute position, select or disable double width, and designate line feed. The printer is preset to work with the default RS-232 interface on Adam balances, but adapts easily to work with other equipment

Print results quickly and easily



Features

- Direct thermal line printing
- Dimensions: 197 x 120 x 102mm (lwxhxh)
- Print size 24 x 24 dot
- Column size 384 dots per line
- 57.5 ± 0.5mm paper width
- Print density: 8 points/mm, 348 points/row
- Print speed of about 50mm per second

Applications

- Printing
- Record-keeping

Accessories

Item number	Description
3126011263	ATP thermal printer paper
3126011281	ATP thermal printer paper (pack of 10)

Key Accessories



ATP Thermal Printer Paper
Print your results and readings with the ATP instead of writing them down. This is replacement paper for the ATP thermal printer. This is one roll of printer paper.



ATP Thermal Printer Paper (pack of 10)
Replacement paper for the ATP thermal printer. This is a 10 pack of printer rolls.

About ADAM:

For 50 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

Supplied by: